



# BESTTRANSFER

## TRANSFER AND REPOSITIONING

# TRANSFER PRODUCTS



## TRANSFERSHEET

Continuous nylon loop fabric with treated inner surface to reduce friction. Glides easily in all directions and significantly reduces force required by caregiver to reposition a patient

part number	size	inches	unit
TS 30110	S	28.5 x 27.5	ea
TS 30120	M	42.5 x 31.5	ea
TS 30130	L	53 x 39.5	ea



## TRANSFERQUILT

Continuous nylon loop fabric with treated inner surface to reduce friction. Glides easily in all directions and significantly reduces force required by caregiver to reposition a patient. Padded for increased patient comfort with same functionality as the transfer sheets.

part number	size	inches	unit
TS 30210	S	23.5 x 23.5	ea
TS 30220	M	23.5 x 55	ea
TS 30230	L	23.5 x 75	ea



## BOOSTERSHEET

Glides easily in up / down direction and significantly reduces force required by caregivers to reposition a patient. Sewn nylon grip handles along both edges assist caregivers.

part number	size	inches	unit
TS 30140	M	29 x 46	ea
TS 30150	L	46 x 46	ea



## DISPOSABLESHEETS

Disposable Transfer sheets have many advantages when it comes to laundry costs and is a great alternative in terms of infectious diseases control.

part number	size	inches	unit
TS 30180	standard	28 x 27.5	50/box



## HANDBELT

Safe lifting belt with non slip inner surface and easy grip handles. Secure buckle system with quick size adjustments. May be worn by patient or caregiver to assist in standing and pivot functions. Superior comfort and leverage versus traditional gait belts.

part number	size	inches	centimeters	unit
TS 30720	M	5 x 40 (adjusts 29 - 45)	12.7 x 101.6 (adjusts 73.7 - 114.3)	ea
TS30730	L	5 x 40 (adjusts 35 x 55)	12.7 x 101.6 (adjusts 88.9 - 139.3)	ea



## HANDIMOVE

Open ended belt with two comfort handles and inner non slip surface. One or two caregivers can apply this by itself or as an aid to sheets/quilts and Handi Lift. Nylon covered with strong poly handles. Inner non slip material covers a foam base.

part number	size	inches	centimeters	unit
TS 30500	standard	9 x 45.5	13 x 115	ea



## HANDISLIDE

Single caregiver aid for quick repositioning. May be used alone or in conjunction with other aids. Bottom surface is low friction while top surface is non slip to adhere to clothing. Strong nylon handle with foam interior.

part number	size	inches	centimeters	unit
TS 30300	standard	9 x 45.5	13 x 115	ea



## HANDILIFT

Mini stretcher for lifting from the floor, lateral transfer or repositioning in bed. Bottom surface is low friction to allow for easy gliding. Four sturdy comfort grip handles allow for use by two or four caregivers as necessary. Maximum weight capacity 300 Lb.

part number	size	inches	centimeters	unit
TS 30400	standard	16.5 x 39.5	42 x 100	ea



## STANDTURN

Two black non slip outer surfaces with slippery inner surfaces. Will easily rotate in all directions when weight applied. User stands in the center and turn in any direction with or without assistance. Designed to be used on beside beds, chairs or outside of cars.

part number	size	inches	centimeters	unit
TS 30820	standard	16	41	ea




## SITTURN

Round double sided disk that spins easily in all directions. Slippery inner surfaces make pivot movements a breeze for people with difficulty turning themselves. Ideal for use on chairs, beds and car seats.

part number	size	inches	centimeters	unit
TS 30810	standard	16	41	ea

## ONEWAYGLIDES SEATING CUSHION

Unique inner surface of these continuous cushioned tubes allows them to glide in only one direction and resist sliding back in the opposite direction. Once placed in position and correctly oriented, the cushion will only move in one direction. When used properly this will help keep a patient correctly seated instead of sliding forward into a slouch position. Promotes good posture and maximizes comfort. Available with a Velour or Non-slip surface.



### NONSLIP

part number	size	inches	centimeters	unit
TS 30610	S	17 x 14.5	43 x 37	ea
TS 30620	M	18 x 18 (best for wheelchairs)	45 x 45	ea
TS 30630	L	31.5 x 23.5	80 x 60	ea



### VELOUR

part number	size	inches	centimeters	unit
TS 30640	S	17 x 14.5	43 x 37	ea
TS 30650	M	18 x 18 (best for wheelchairs)	45 x 45	ea
TS 30660	L	31.5 x 23.5	80 x 60	ea