BESTMOVETM


STA400 | STANDING TRANSFER AID

The BestMove[™] STA400 is your standing transfer solution. Transfers are quick and require minimal assistance. Fold-up seats and an easy grip crossbar allows the user to actively participate in a transfer process. Tapered legs allow unit to squeeze in between tight spaces.




FEATURES

- 400 lb capacity
- Shin brace for secure transports
- Support strap available attachment hooks for optional side support and back strap
- Affordable alternative to electric stand assists
- First choice for any fall prevention program



SPECIFICATIONS FOR STA400	
maximum capacity	400Lb / 180Kg
min. base inner width	13.7"
max. base outer width	20''
fixed base outer back	25"
width	
seat width	20''
overall length	31"
seat height	26.3"
overall height	43''

^{*} SPECIFICATIONS PROVIDED ARE A GENERAL GUIDELINE AND ARE SUBJECT TO CHANGE WITHOUT NOTICE



20"