BESTMOVETM

STA450 | STANDING TRANSFER


The BestMove[™] is your best choice for fall prevention and is also a vital part of reducing lifting injuries while at the same time promoting mobility for the patient. Offers ample weight capacity at 450 lb max. Includes an adjustable base to provide clearance around common obstacles such as recliners allowing access to patients.



| FEATURES

- 450 lb. Capacity
- Dual shin pads for extra support
- Optional back support strap available
- Simple solution to compliment your electric stand assist lifters

DIMENSIONS



SPECIFICATIONS FOR STA450	
maximum capacity	450 lb. / 205 kg
minimum base width	26.8"
open base width	37''
overall height	43.3"
minimum seat height	26.7"
seat width	20''
overall length	35.4"
shin pads	yes
optional support strap	yes
optional seat locks	yes

^{*} SPECIFICATIONS PROVIDED AS GENERAL GUIDELINE AND ARE SUBJECT TO CHANGE WITHOUT NOTICE

